

YieldAce Hotstart DNA Polymerase

INSTRUCTION MANUAL

Catalog #600334 (1000 U) and #600336 (3000 U)
Revision A.01

For In Vitro Use Only
600330-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

YieldAce Hotstart DNA Polymerase

CONTENTS

Materials Provided	1
Storage Conditions	1
Notices to Purchaser	1
Introduction	2
Critical Optimization Parameters	3
DNA Template	3
Enzyme Concentration	3
Deoxynucleotide Concentrations.....	4
Primer Design and Concentration	4
Cycling Parameters.....	4
Protocol	5
Troubleshooting	7
References	8
Endnotes	8
MSDS Information	8

YieldAce Hotstart DNA Polymerase

MATERIALS PROVIDED

Materials provided	Catalog #600334 (1000 U)	Catalog #600336 (3000 U)
YieldAce hotstart DNA polymerase (5 U/ μ l)	200 μ l	3 \times 200 μ l
10 \times YieldAce reaction buffer	4 \times 1 ml	12 \times 1 ml

STORAGE CONDITIONS

All Components: -20°C

NOTICES TO PURCHASER

Notice to Purchaser: Limited License

Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights (such as 5' Nuclease Process patent rights) are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

INTRODUCTION

YieldAce Hotstart DNA Polymerase* is a unique enzyme formulation specifically developed for maximizing PCR product yields while providing hotstart capabilities.¹ YieldAce Hotstart DNA polymerase is formulated with heat labile monoclonal antibodies that, at room temperature, effectively neutralize DNA polymerase activity. Full enzyme activity is regained upon denaturation of the antibody during the initial denaturation step. The enzyme and PCR reaction buffer have been optimized to provide superior performance in a variety of applications routinely carried out with *Taq* DNA polymerase. Since YieldAce DNA polymerase differs significantly from *Taq* DNA polymerase, it is important to follow the recommended protocol for optimal performance. Compared to *Taq* DNA polymerase, YieldAce DNA polymerase synthesizes higher yields of product from a variety of DNA templates, including plasmid DNA, genomic DNA, and cDNA, and it readily amplifies templates directly from bacterial cultures. Additionally, reactions performed with YieldAce DNA polymerase show less variation in product yield, with respect to amplicon size and DNA template concentration, than identical reactions performed with *Taq* DNA polymerase.

A key component of YieldAce hotstart DNA polymerase is the ArchaeMaxx polymerase-enhancing factor. The ArchaeMaxx factor eliminates a PCR inhibitor and promotes shorter extension times, higher yield, and greater target length capabilities. The ArchaeMaxx factor improves the yield of products by overcoming dUTP poisoning, which is caused by dUTP accumulation during PCR through dCTP deamination.² Once incorporated, dU-containing DNA inhibits *Pfu* and most archaeal proofreading DNA polymerases, such as Vent® and Deep Vent® DNA polymerases, limiting their efficiency.² The ArchaeMaxx factor functions as a dUTPase, converting poisonous dUTP to harmless dUMP and inorganic pyrophosphate, resulting in improved overall PCR performance.

* U.S. Patent Nos. 6,734,293; 6,489,150, 6,444,428, 6,183,997, 5,948,663, 5,866,395, 5,556,772, 5,545,552 and patents pending.

CRITICAL OPTIMIZATION PARAMETERS

All PCR amplification reactions require optimization to achieve the highest product yield and specificity. Critical optimization parameters for successful PCR using YieldAce hotstart DNA polymerase are outlined in Table I and discussed in the following section. These include the quality and concentration of DNA template, the appropriate concentration of dNTPs, sufficient enzyme concentration, optimal cycling parameters, and appropriate PCR primer design. The magnesium chloride supplied in the 10× reaction buffer is optimized for YieldAce hotstart DNA polymerase and most targets. Further optimization is rarely necessary.

TABLE I
Optimization Parameters and Suggested Reaction Conditions

Parameter	Suggested Condition (50 µl reaction volume)
Input template	100–250 ng genomic DNA 1–30 ng vector DNA 0.5–2 µl bacterial culture 1–4 µl cDNA synthesis reaction
YieldAce hotstart DNA polymerase	2.5 U
dNTP concentration	200 µM each dNTP (0.8 mM total)
Primers (each)	100 ng (0.25 µM)
Extension time	1 min per kb
Denaturing temperature	95°C
Extension temperature	72°C

DNA Template

Template Quality

Potential shearing of the genomic DNA template is minimized by the use of wide-bore tips for pipetting or mixing of the template. Freezing of high molecular weight templates should be avoided; storage at 4°C is recommended. Bacterial cultures may be directly amplified from glycerol stocks or fresh overnight cultures.

Template Concentration

When amplifying genomic DNA templates, use 100–250 ng of template for reaction volumes of 50 µl. To amplify low-complexity targets (for example, lambda or plasmid DNA), use 1–30 ng of template or 0.5–2 µl of a bacterial culture. To amplify cDNA, 1–4 µl of a standard cDNA reaction is recommended.

Enzyme Concentration

Robust product yield requires an adequate DNA polymerase concentration. The use of 2.5 U/50-µl reaction consistently generates high yield from most templates. Higher concentrations of enzyme in the reaction mixture can be inhibitory.

Deoxynucleotide Concentrations

Amplification efficiencies are influenced by deoxynucleotide (dNTP) concentrations. Insufficient concentrations of dNTPs may result in lower yields. We recommend using each dNTP at a final concentration of 200 μM .

Primer Design and Concentration

Primers should be ≥ 23 bp in length with a balanced $T_m \geq 60^\circ\text{C}$. The resulting high annealing temperature promotes specificity and discourages secondary structure formation. Further, primer sequences should be analyzed for potential duplex and hairpin formation as well as false priming sites in order to obtain the highest yield of specific PCR products.

We recommend using ~ 0.25 μM final concentration of each primer. When using 25-mer oligonucleotide primers in a 50- μl reaction volume, this is equivalent to ~ 100 ng of each primer.

Cycling Parameters

As with all PCR reactions, cycling parameters are critical for successful amplification and may require further optimization.

Extension Time

Maintain an extension time of 1.0 minute (for targets ≤ 1 kb) or 1.0 minute/kb of template (>1 kb targets) for general applications; longer extension times may produce higher yields.

Extension Temperature

Extension temperatures also have a critical effect on amplicon yield. An extension temperature of 72°C should be used with all templates.

Denaturation Temperature

High denaturation temperatures damage DNA templates, so the denaturation temperature should be as low as possible. A denaturation temperature of 95°C works well for most targets.

PROTOCOL

1. Prepare a reaction mixture for the appropriate number of samples to be amplified. The following table provides an example of reaction mixtures for the amplification of genomic and plasmid DNA targets. The recipe listed in the table is for one reaction and can be scaled for multiple samples. Add the components *in order* and mix gently.

Component	Quantity per reaction	
	Genomic targets	Plasmid targets
Distilled water	X μ l to final volume of 50.0 μ l	X μ l to final volume of 50.0 μ l
10 \times YieldAce reaction buffer	5.0 μ l	5.0 μ l
dNTP mix (25 mM of each dNTP)	0.4 μ l	0.4 μ l
DNA template: Genomic DNA Plasmid DNA	100–250 ng* —	— 1–30 ng* or 0.5–2 μ l bacterial culture
Primer #1	100 ng	100 ng
Primer #2	100 ng	100 ng
YieldAce hotstart DNA polymerase (5 U/ μ l)	0.5 μ l	0.5 μ l
Total reaction volume	50.0 μ l	50.0 μ l

* Excess template can be inhibitory to the PCR reaction. Do not exceed the recommended quantities.

2. Before thermal cycling, aliquot 50 μ l of the master mixture into sterile thin-walled PCR tubes.
3. If the temperature cycler is not equipped with a heated cover, overlay each reaction with ~50 μ l of DNase-, RNase-, and protease-free mineral oil (available from Sigma Chemical Company, St. Louis, Missouri).

4. Perform PCR using optimized cycling conditions. Suggested cycling parameters are given below for (A) single-block temperature cyclers and (B) RoboCycler temperature cyclers:

(A)

Single-Block Temperature Cyclers

Segment	Number of cycles	Temperature	Duration
1	1	92°C	2 minutes
2	10	95°C	20 seconds
		Primer $T_m - 5^\circ\text{C}^a$	20 seconds
		72°C	60 seconds/kb of PCR target ^b
3	20	95°C	20 seconds
		Primer $T_m - 8^\circ\text{C}^a$	20 seconds
		72°C	60 seconds/kb of PCR target plus 10 seconds/cycle ^b
4	1	72°C	7 minutes

^a The annealing temperature may be lowered or raised if necessary to obtain optimal results.

^b The minimum extension time is 1 minute.

(B)

RoboCycler Temperature Cyclers

Segment	Number of cycles	Temperature	Duration
1	1	92°C	2 minutes
2	10	95°C	50 seconds
		Primer $T_m - 5^\circ\text{C}^a$	50 seconds
		72°C	90 sec./kb ^b
3	20–25	95°C	50 seconds
		Primer $T_m - 8^\circ\text{C}^a$	50 seconds
		72°C	90 sec./kb ^b
4	1	72°C ^b	7 minutes

^a The annealing temperature may be lowered or raised if necessary to obtain optimal results.

5. Analyze the PCR amplification products by electrophoresis using an appropriate percentage acrylamide or agarose gel.

TROUBLESHOOTING

Observations	Suggestions
No PCR product or lower yield than expected	Verify the concentration of DNA template. Titrate the amount of full-length intact DNA template and/or increase the number of cycles up to a maximum of 40 cycles
	Use intact and highly purified DNA templates
	Store the template at 4°C; do not freeze the template
	Lower the annealing temperature in 5°C increments
	Allow at least 60 seconds of extension time for each kilobase to be amplified (90 seconds of extension time per kilobase may also be helpful for difficult templates)
	Denaturation times of 50 seconds for RoboCycler temperature cyclers or 20 seconds for single-block temperature cyclers at 95°C are usually sufficient, longer denaturation times or higher denaturation temperatures may damage the DNA template; use the shortest denaturation time compatible with successful PCR on the thermal cycler
	Primer pairs exhibiting matched primer melting temperatures (T_m) and complete complementarity to the template are recommended
	Analyze the primer sequences to ensure that duplexes or hairpins do not form
	Gel-purified or HPLC-purified primers ≥ 23 nucleotides in length are desired for successful PCR
	Adjust the ratio of primer versus template to optimize yield of the desired product
	Use thin-wall PCR tubes for Stratagene thermal cyclers. These PCR tubes are optimized to ensure ideal contact with the multiblock design to permit more efficient heat transfer and to maximize thermal-cycling performance
	Use DMSO in the PCR mixture; titrate the DMSO concentration in 1% increments
Artifactual PCR smears	Decrease the amount of YieldAce Hotstart polymerase
	Reduce the extension time
	Optimize the cycling parameters specifically for the primer-template set and the thermal cycler used
Multiple bands	Increase the annealing temperature in 3°C increments in segment 2
	Use Perfect Match PCR enhancer to improve PCR product specificity
	Use DMSO in the PCR mixture; titrate the DMSO concentration in 1% increments
	Verify that the primers hybridize only to the desired sequences on the template

REFERENCES

1. Borns, M. and Hogrefe, H. H. (2000) *Strategies* 13(4):125–127.
2. Hogrefe, H. H., Hansen, C. J., Scott, B. R. and Nielson, K. B. (2002) *Proc Natl Acad Sci U S A* 99(2):596–601.
3. Borns, M., Scott, B. and Hogrefe, H. H. (2001) *Strategies* 14(1):5–8.
4. Eads, J., Bai, F., Allen, R. and Hogrefe, H. H. (2000) *Strategies* 13(4):138–140.

ENDNOTES

Deep Vent® and Vent® are registered trademarks of New England Biolabs, Inc.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.